



Stent Graft System Overview

STEP 1 Select ALTO Aortic Body

Stent Graft Diameter (mm)	Aortic ID* (mm)
20	16 - 17
23	18 - 20
26	21 - 23
29	24 - 26
34	27 - 30



*At the intended proximal sealing ring location.

STEP 2 Select Ovation iX Iliac Limbs

Labeled Diameter (mm)	Native Iliac Vessel Id Range (mm)
10	8 - 9
12	10 - 11
14	12 - 13
16	14 - 15
18	16 - 17
22	18 - 20
28	21 - 25

STEP 3 Select Ovation iX Extensions

Alto[®] Abdominal Stent Graft System



Stent Graft System Specifications

ALTO AORTIC BODY	Catalog Number	Delivery System Inner Profile (F)	Delivery System Outer Profile (F)	Stent Graft Proximal Diameter (mm)		Stent Graft Length (mm)	
	TV-AB2080-N	13	15	20		80	
TV-AB2380-N	13	15	23		80		
TV-AB2680-N	13	15	26		80		
TV-AB2980-N	13	15	29		80		
TV-AB3480-N	13	15	34		80		

OVATION IX ILIAC LIMB	Catalog Number	Delivery System Inner Profile (F)	Delivery System Outer Profile (F)	Stent Graft Proximal Diameter (mm)	Stent Graft Distal Diameter (mm)	Stent Graft Length (mm)	Stent Graft Length with Aortic Body (mm)
	TV-IL141080-J	10	12	14	10	80	125
TV-IL1410100-J	10	12	14	10	100	145	
TV-IL1410120-J	10	12	14	10	120	165	
TV-IL1410140-J	10	12	14	10	140	185	
TV-IL1410160-J	10	12	14	10	160	205	
TV-IL141280-J	10	12	14	12	80	125	
TV-IL1412100-J	10	12	14	12	100	145	
TV-IL1412120-J	10	12	14	12	120	165	
TV-IL1412140-J	10	12	14	12	140	185	
TV-IL1412160-J	10	12	14	12	160	205	
TV-IL141480-J	10	12	14	14	80	125	
TV-IL1414100-J	10	12	14	14	100	145	
TV-IL1414120-J	10	12	14	14	120	165	
TV-IL1414140-J	10	12	14	14	140	185	
TV-IL1414160-J	10	12	14	14	160	205	
TV-IL141680-J	11	13	14	16	80	125	
TV-IL1416100-J	11	13	14	16	100	145	
TV-IL1416120-J	11	13	14	16	120	165	
TV-IL1416140-J	11	13	14	16	140	185	
TV-IL1416160-J	11	13	14	16	160	205	
TV-IL141880-J	11	13	14	18	80	125	
TV-IL1418100-J	11	13	14	18	100	145	
TV-IL1418120-J	11	13	14	18	120	165	
TV-IL1418140-J	11	13	14	18	140	185	
TV-IL1418160-J	11	13	14	18	160	205	
TV-IL142280-J	12	14	14	22	80	125	
TV-IL1422100-J	12	14	14	22	100	145	
TV-IL1422120-J	12	14	14	22	120	165	
TV-IL1422140-J	12	14	14	22	140	185	
TV-IL1422160-J	12	14	14	22	160	205	
TV-IL142880-J	13	15	14	28	80	125	
TV-IL1428100-J	13	15	14	28	100	145	
TV-IL1428120-J	13	15	14	28	120	165	
TV-IL1428140-J	13	15	14	28	140	185	
TV-IL1428160-J	13	15	14	28	160	205	

OVATION IX ILIAC EXTENSION	Catalog Number	Delivery System Inner Profile (F)	Delivery System Outer Profile (F)	Stent Graft Proximal Diameter (mm)	Stent Graft Distal Diameter (mm)	Stent Graft Length (mm)
	TV-EX101045-J	10	12	10	10	45
TV-EX121245-J	10	12	12	12	45	
TV-EX141445-J	10	12	14	14	45	
TV-EX161645-J	11	13	16	16	45	
TV-EX181845-J	11	13	18	18	45	
TV-EX222245-J	12	14	22	22	45	
TV-EX282845-J	13	15	28	28	45	

CUSTOMSEAL™ POLYMER FILL KIT

Catalog Number
TV-CS14-G



AUTOINJECTOR 2

Catalog Number
TV-AI01-N



The ALTO® Abdominal Stent Graft System and associated components are not available in all countries or regions. Please contact your Endologix representative for details regarding product availability. Prior to use, refer to Instructions for Use for more information concerning Indications, Contraindications, Specific Anatomic Considerations, Warnings, Precautions, and Adverse Events.

Rx only.

©2021 Endologix LLC. All rights reserved. MM2148 Rev 05

